

### DURABLE and SUSTAINABLE CENTRIFUGES

#### MAIN ADVANTAGES:

- > Solid metal housing
- > Lid Locking System – SAFE TO USE
- > Variable Speed with Closed Loop Regulation
- > Electronic Timer with HOLD function
- > Microprocessor Controlled Functions of Cover, Time and Speed



Clinical centrifuge Centric 90 is appropriate for veterinary practice and other operations in small laboratories. It is designed for optimal and efficient separation of blood and other materials with different densities. Irrespective of simple construction Centric 90 has many advantages and carries a CE mark.

TECHNICAL DATA	CENTRIC 90
Code for 230V 115V	456.900 1456.900
Power supply	230V / 50-60Hz 115V / 50-60Hz
Power requirement	100W
Fuses	2x1A T - 230V 2x2A T - 115V
Rotor	Angle rotor RA6 with 12 or 15ml tubes
Max. rotation speed	500 to 3.500 rpm, in steps of 100 rpm
Max. centrifugal force	1300 x g / 12ml 1600 x g / 15ml
Time	0 - 30 min, HOLD function
Permitted ambient temperature	0°C to + 40°C
Dimensions (H x W x D)	340 x 290 x 225 mm
Weight	11kg

*We reserve the right to alter specification details etc. without prior notice or liability!*